
Element Biosciences

OFFICIAL INDIA CHANNEL PARTNER
Genique Lifesciences™

Cytos™ Cell Screening

See the biology driving every screen




Morphology

Provides robust segmentation and links perturbations to cell state




Transcriptome

Resolve gene program and pathway changes




Profile morphology, transcriptome, protein, and guide RNAs from the same cell in one automated run with Cytos on AVITI24™.



Guide RNA

Directly sequencing 10s-10ks of guide RNAs for genome-wide or focused screens.



Protein

Visualize protein expression and localization

CRISPR screening

Turn every genetic perturbation into a multidimensional map. Sequence guide RNA with morphology, protein, and 3′ transcriptome for a complete view of biology.

Readouts: guide RNA · morphology · up to 20-plex protein · optional 3′ transcriptome

Drug screening

Connect a compound's effect to its mechanism with same cell multimodal readouts in every run. Understand how a drug changes cell state, protein localization or expression, and molecular pathways.

Readouts: 3′ transcriptome · morphology · up to 20-plex protein


Pre-order Cytos now

Be among the first labs running Cytos.

For Research Use Only. Not for use in diagnostic procedures. | August 2026


Streamlined workflow from cells to answers

Cytos uses DISS, Direct In Sample Sequencing, to capture every readout from the same fixed, intact cell, with automated analysis.




Sample prep

Culture cells directly on Cytos slides, apply treatments, fix, and load. No library prep.



Imaging & sequencing

Morphology, protein, transcriptome, and guide assignment captured on AVITI24.



Secondary analysis

ElemBio Catalyst™ handles transcriptome processing immediately following your run.



Visualization and tertiary analysis

Bring every readout into one view for QC and downstream multimodal analysis. Data compatible with your existing multiomic pipelines.

What makes Cytos different?

Same-cell multiomics

No separate assays to reconcile.

One automated instrument

All your readouts, start to finish, on AVITI24.

Easy customization

Choose custom morphology and protein markers that fit your screen.

No library prep

DISS sequences directly in intact cells with no separate NGS run.

Speed to insight

Less than 1 hour of offboard prep and run times as fast as 1 day.

End-to-end workflow

From sample to visualized result.

Information in this document is provided for research use only and is subject to change without notice.
© 2026 Element Biosciences, Inc. All rights reserved. AVITI24, ElemBio Catalyst, and Cytos are trademarks and/or registered trademarks of Element Biosciences, Inc. in the United States and/or other countries. All other trademarks, service marks, and logos are the property of their respective owners. Use of third-party marks is for identification purposes only and does not imply endorsement, sponsorship, or affiliation.

Visit elementbiosciences.com for more information.

Pre-order Cytos now

Be among the first labs running Cytos.

OFFICIAL INDIA CHANNEL PARTNER

Genique Lifesciences™

contact@genique.co | +91 124 426 3051

Element Biosciences